

## 2. Experimental procedure

### OBLIGATORY:

Please watch on YouTube the clip "Automatic Pipet | MIT Digital Lab Techniques Manual" explaining how to use automatic pipets. You can find it under this link:

<https://youtu.be/AcNtVgOp0bl>

### Material:

- Competent cells of *E. coli* M15 (pREP4)
- Plasmid pQE32-Eos
- Sterile LB (culture medium for bacteria)
- Sterile agar
- Antibiotic solution (ampicillin 100 mg/ml, kanamycin 50 mg/ml)
- IPTG (1 M)
- Ni<sup>2+</sup>-charged resin
- Chromatography column
- Buffers for protein purification:
  - Lysis buffer (20 mM Tris-HCl pH 7.5, 300 mM NaCl)
  - Elution buffer (20 mM Tris-HCl pH 7.5, 300 mM NaCl, 300 mM Imidazole)

Wear gloves and a lab coat!

### 2.1 Transformation of competent *E. coli* M15 (pREP4) cells with plasmid pQE32-Eos

#### Day 1:

- Thaw two aliquots (one sample, one negative control) of competent *E. coli* M15 (pREP4) cells on ice for 5 min.
- Add 1 µl of the plasmid to the sample cells and mix by pipetting up and down.
  - No plasmid is added to the negative control.
- Incubate both aliquots on ice for 5 min.
- Heat shock both aliquots for 45 sec at 42 °C.
- Incubate both aliquots on ice for 5 min.
- Add 1 ml of sterile LB to each aliquot.
- Incubate for 1 h shaking at 37 °C.
- In the meantime, prepare two agar plates and label them: Sample and negative control.
- Plate 100 µl of the sample and the negative control on one agar plate each.
- Incubate the plates over night at 37 °C.

#### Preparation of two agar plates:

- Melt 25 ml 2x agar in the microwave and mix with 25 ml 2x LB.
- Cool the mixture down to about 40 °C.
- Add 50 µl of ampicillin and 25 µl of kanamycin.
- Mix by pipetting up and down.
- Pour the medium in 2 different agar plates.

## 2.2 Bacterial growth and Eos overexpression

### Day 2 morning:

5 ml preculture:

- Mix 5 ml sterile LB with 5  $\mu$ l ampicillin and 2.5  $\mu$ l kanamycin.
- Inoculate the medium with bacteria grown on the agar plate: Scratch the plate using a sterile pipet tip and drop the tip into the test tube.
- Let the bacteria grow at 37 °C and 200 rpm for 6 h.

### Day 2 afternoon:

100 ml culture:

- Prepare 100 ml sterile LB with 100  $\mu$ l ampicillin and 50  $\mu$ l kanamycin.
- Inoculate the medium with the bacteria from the preculture ( $OD_{600}=0.06$ ).
- Cultivate the cells at 37 °C and 200 rpm shaking and measure the  $OD_{600}$  every 20 min (show the growth curve in your report).
- When  $OD_{600}$  is at 0.5 - 0.7, add 100  $\mu$ l IPTG to the culture.
- Grow over night at 25 °C and 200 rpm shaking.

## 2.3 Cell lysis and protein purification

### Day 3 morning:

Harvesting of the cells and cell lysis:

- Harvest the bacteria culture by centrifugation (10 min, 3000-6000x g, 4 °C)
- Discard the supernatant and put the pellet on ice.
- Add 30 ml lysis buffer to the pellet and resuspend the pellet in the buffer.
- Lyse the cells by sonication on ice:
  - Use the Branson Sonifier<sup>®</sup> in lab N25/3204: 1 min effective sonication time with pulses of 15 sec at 30 % power (amplitude) with resting on ice for 30 sec between the pulses. *no We also covered it to prevent photobleaching.*
- Centrifuge the sample (20 min, 20000x g, 4 °C). The supernatant (cytosol) should contain the protein. *no lunch break*

Protein purification (Nickel affinity chromatography):

- Wash 0.5 ml mixed  $Ni^{2+}$ -charged resin per group once with water and once with lysis buffer (centrifugation for 2 min at 1000x g).
- Add the cytosol to the resin and incubate it in a rotating wheel (10 min, 4 °C).
- Centrifuge (2 min, 1000x g) and carefully remove the supernatant.
- Resuspend the resin in 0.5 ml lysis buffer, transfer it to a spin column by pipetting and then centrifuge (1 min, 1000x g).
- Wash the resin with 0.8 ml lysis buffer and centrifuge (1 min, 1000x g).
- Remove residual wash buffer by centrifuging again (1 min, 1000x g).

- Transfer the spin column to a fresh tube, add 0.8 ml elution buffer to the resin, incubate for 2 min and then centrifuge (1 min, 1000x g).

## 2.4 Protein quantification and verification

Day 3 afternoon:

Computed  
Recording of the absorption spectrum and calculation of the protein concentration:

$$\epsilon_{280} = 27,452.5$$

$$\epsilon_{506} = \pm 62.5$$

$$72,000$$

$$M^{-1}cm^{-1}$$

Prior to the experiment: Look up the extinction coefficient of Eos at 280 nm and 506 nm.

- Measure the absorption spectrum of the protein solution using a UV spectrometer (Cary Varian 50) in the range 250 nm-600 nm,  $OD_{506} < 1$ . Use the 'Advanced Settings': Scan rate: medium. First, scan the elution buffer as a blank (baseline).
- Calculate the protein concentration using the molecular extinction coefficient of Eos (use Beer-Lambert law).
- Dilute the protein such that the final concentration is 1  $\mu M$  in 200  $\mu l$  final volume.

Day 4:

Photoconversion of Eos:

**CAUTION:** The laser with a power of up to 100 mW can irrevocably damage the human eye! Therefore, it must be operated only after the lid of the sample chamber was closed!

In order to verify that the product is in fact Eos, an experiment is performed, which shows that the green-emitting protein is converted to a red-emitting protein upon illumination with UV light.

- Fill the above-prepared solution ( $\approx 1 \mu M$  Eos) into a special black fluorescence cuvette (1 cm path length, 9 mm beam height, open side window).
- Place the cuvette into the photoconversion setup for monitoring fluorescence from both the green and the red form of Eos.
- Start OceanView software and open project 'Fast\_Photoconv', which contains the following settings:
  - 'integration time': 50 ms (reduce, if saturation ( $>60,000$ ) occurs)
  - 'scans to average': 20 (update time is then 1 second)
  - 'boxcar width': 5 (corresponds to  $\approx 1$  nm spectral resolution)
- Configure graph saving (file symbol with tool on the right). Most settings are already done in the project 'Fast\_Photoconv' but the bold items have to be set manually!
  - **Set Target directory: group#** (directory has to be generated before)
  - File format: ASCII (plain)
  - File suffix: file counter
  - Save: every scan
  - Stop after this many scans: 60
  - Apply. Exit.

No other Group

Soban - 85,000

78,000 43,000

- Rescale the graph to the range 400-700 nm (magnifier symbol with the '2' inside)
- **Push the rightmost button 'view trendline time trail'**
- When the cuvette is in place:
  - Start data acquisition by the 'start all file writers' symbol and close the lid of the sample chamber, which closes a contact to **switch on the laser** (4.5 V corresponding to  $\approx 33$  mW of laser power)
  - 'File saving in progress' will appear at the lower left
  - (Note: most of the settings can be changed during acquisition, even the stopping conditions!)
  - Keep the box closed while the laser is running**
- When the measurement is finished:
  - Move the series of 60 spectra to a USB stick (2-column text-files)
  - Produce at home a graph showing a few spectra ( $\sim 5-10$ ) at different stages with corresponding times indicated.
  - Produce a kinetics graph (intensity over time) for the peak intensity of the green and the red species.

## LITERATURE:

- [1] Alberts, B. *et al.*: "Molecular Biology of the Cell", Garland Science, 5<sup>th</sup> edition
- [2] Lakowicz, J.R.: "Principles of Fluorescence Spectroscopy", Kluwer Academic/Plenum Publishers
- [3] Sambrook, J.: "Molecular Cloning", Cold Spring Harbour Laboratory Press, 2000.
- [4] Wiedemann, J. *et al.*: "EosFP, a fluorescent marker protein with UV-inducible green-to-red fluorescence conversion", PNAS, 2004
- [5] Nienhaus, G.U. *et al.*: "Photoconvertible Fluorescent Protein EosFP", Biophysical Properties and Cell Biology Applications. Photochemistry and Photobiology, 2006.



1°	0.107	00
2°	0.130	
3°	0.205	
4°	0. <del>244</del> 244	
5°	0.448	
6°	0.675	
7°		

We use position n° 2